

# Tush Push Toilet Seat Lift Installation Guide

Note: You will need a 5/8" and a 9/16" wrench for installation

Note: The toilet in the diagrams (except for the first one) does not have a tank installed.


1. Remove the existing seat from your toilet.
2. Install the main toilet lift frame on the toilet bowl. The mounting screws will fit through the bowl. Attach the plastic bolts to the bottom side of the screws.



3. (3a) Lift up on the top portion of the frame from the back, it should stay up for the placement of the motor mount and motor. (3b) Place the motor mount over the bolt.



4. Push the motor mount down behind the frame bracket and it is ready to mount the motor.



# Tush Push Toilet Seat Lift Installation Guide

5. Place the bottom motor mount slot over the lower motor mount bracket and push in the pin. Hold the top shaft of the motor and lower the top portion of the main frame. Align the motor slot with the top frame bracket and push in the pin. This completes the installation.



6. Place the two button hand set on the arm and plug in the unit.
7. Adjust arms to desired position and tighten knobs.

## Cleaning Instructions

1. With the handset, raise the toilet lift to its "up" position.
2. Lift toilet seat to allow access to removable shield.
3. Grab the front of the shield. (See Figure 1).
4. Pull upward and then forward until the shield slides off of support rollers (See Figure 2).
5. Clean and reinstall shield. Use a cloth with a mild cleaner. **DO NOT USE TOILET BOWL CLEANER OR A SPRAY BOTTLE.**



### Replacement Shield

If your removable shield is showing signs of wear or gets damaged, we offer replacement shields at a reasonable price. Log-on to our website [www.phillipsliftsystems.com](http://www.phillipsliftsystems.com) and order online, or call us at: 888-429-5438 and we'll ship a new replacement shield to you that same day.