

B A R I A T R I C P O W E R C H A I R



# Roll

Seat Colors  
Charcoal/Gray  
Coffee/Sand

Body Colors  
Candy Apple Red  
Blue

*Alanté II*  
BY GOLDEN

Heavy duty rear wheel drive power chair offering an oversized seat, rugged motors, and a powerful controller.



YOUR LIFE IN MOTION

**GOLDEN**  
technologies

## Alanté GP-202 Specifications

<b>Weight:</b>	
Base Frame	91 lbs
Seat w/ arms/controller	48 lbs
<b>Overall Dimensions:</b>	
Length	40"
Width	27"
Seat Height	20" to 23" – adjustable
<b>Top Speed:</b>	5 mph
<b>Operating Range:</b>	12 miles*
<b>Load Capacity:</b>	450 lbs
<b>Turning Radius:</b>	30.5"
<b>Wheel Sizes:</b>	
Drive Wheels	10"
Casters	8"
Anti-tippers	6"
<b>Seat Dimensions:</b>	
Headrest Height	adjustable
Depth	20"
Width	20"
<b>Footrest:</b>	
Height Adjustable	Yes
Angle Adjustable	Yes
<b>Controller:</b>	Dynamic Shark 75 amp, programmable
<b>Charger:</b>	Off Board
<b>Color:</b>	Blue, Candy Apple Red
<b>Seat Color</b> (standard & high back):	Coffee/Sand, Charcoal/Gray
<b>Battery Type:</b>	2 (two) 22NF 12V



- Four point seat system
- Multiple seat height positioning
- Rear power mount for easy connections

## Options

- Elevating Leg Rests
- Specialty Leg Supports
- Custom Seating
- Swing Away Footrests

## Features

- Adjustable 4" anti-tip wheels
- Powerful battery charger
- Adjustable foot plate
- Two brilliant colors
- 450 lb. weight capacity
- Adjustable head rest
- Adjustable armrests
- Adjustable seat (pitch/height)
- Strong, rugged motors

\* Battery range will vary due to rider weight, drive surface, drive slope, and tire pressure.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.