


Sunrunner™ 4

(888B-4)

The **Sunrunner™ 4** comes equipped with a deluxe captain seat, full tiller shroud, wicker basket and mirror. With four wheels, this scooter can move and turn securely. The **Sunrunner™ 4** is perfect for those who want stability and luxury in their everyday life at a value price.



Available Colors



- **300 lbs. weight capacity**
- **Sliding and swivel seat**
- **Mirror and adjustable tiller**
- **Deluxe captain seat and full tiller shroud**
- **Headlight package**
- **Front chrome bumper**

Overall Dimensions (L' x W x H)	in	49 x 22 x 38
Number, Size of Tires (front; rear)	Front (rear)	2, 10" (2, 10")
Suggested User Weight (on a level road)	lb	Max. 300
Weight of Heaviest Part (w/o baskets and seat)	lb	N/A
Total Weight (w/ batteries, basket, seat)	lb	180
Each Battery Module Weight	lb	24
Battery Module Capacity	_V_Ah x Pcs	12V33-38Ah x 2
Charger		Off Board (3A)
Output Power of Motor	hp	1
Maximum Speed ²	mph	5
Max. Climbing Angle/Safe Climbing Angle (based on user weight of 90kg/200lb.)	% (°)	14 (8)
Range (per charge) (After the battery and mechanical moving parts fully break in)	mile	25
Turning Radius	in	50
Ground Clearance	in	3
Handle Bar	Type	T-Bar

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

b) Ground surface

c) Battery capacity and conditions

d) Type of charger

e) Ambient temperature

f) The way of driving

g) If the battery and mechanical

moving parts fully break in

h) Etc.

SHOPRIDER

Toll-Free: 800.743.0772

www.shoprider.com