


Sunrunner™ 3

(888B-3)

One of Shoprider's most popular models is the **Sunrunner™ 3**. This scooter is full of features including slide and swivel seating, front headlight package and much more at a great value. The **Sunrunner™ 3** is perfect for maneuvering indoors and outdoors.



Available Colors



- 300 lbs. weight capacity
- Sliding and swivel seat
- Adjustable tiller
- Deluxe captain seat and full tiller shroud
- Headlight package
- Seat height adjustable
- Basket included
- Rear light package w/ blinkers/ hazards upgrade avail.
- Rear view mirror included

Overall Dimensions (L x W x H)	in	45 x 24 x 38
Number, Size of Tires (front; rear)	Front (rear)	1, 9" (2, 10")
Suggested User Weight (on a level road)	lb	Max. 300
Weight of Heaviest Part (w/o baskets and seat)	lb	N/A
Total Weight (w/ batteries, basket, seat)	lb	165
Each Battery Module Weight	lb	24
Battery Module Capacity	_V_Ah x Pcs	12V33-38Ah x 2
Charger		Off Board (3A)
Output Power of Motor	hp	1
Maximum Speed ²	mph	5
Max. Climbing Angle/Safe Climbing Angle (based on user weight of 90kg/200lb.)	% (°)	14 (8)
Range (per charge) (After the battery and mechanical moving parts fully break in)	mile	25
Turning Radius	in	37
Ground Clearance	in	3
Handle Bar	Type	T-Bar

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

b) Ground surface

c) Battery capacity and conditions

d) Type of charger

e) Ambient temperature

f) The way of driving

g) If the battery and mechanical

moving parts fully break in

h) Etc.

SHOPRIDER

Toll-Free: 800.743.0772

www.shoprider.com