

SPRINTER™ JUMBO XL 889X

Shoprider is proud to introduce the revolutionary **Sprinter Jumbo XL**. With all new 14" five-spoke sport wheels, 23" black leatherette pillow top seat and four-wheel suspension, the **Sprinter Jumbo XL** is truly ahead of its time. This bold, daring and sporty new scooter is perfect for those who are looking for comfort and style with added weight capacity.

Available Colors



Burgundy



- **500 lb. weight capacity**
- **Slide and swivel seat**
- **23" luxurious leatherette executive seat**
- **Full light package**
- **Four wheel suspension**
- **All new 14" five-spoke sport wheels**

Overall Dimensions (L1 x W x H)	in	56 x 25 x 37
Number, Size of Tires (front, rear)	Front (Rear)	2, 14" (2, 14")
Suggested User Weight (on a level road)	lb (kg)	Max. 500 (227)
Weight of Heaviest Part (w/o baskets and seat)	lb (kg)	107 (49)
Total Weight (w/ batteries, basket, seat)	lb (kg)	208 (70)
Each Battery Module Weight	lb (kg)	32 (15)
Battery Module Capacity	_V_Ah x Pcs	12V33-75Ah x 2
Charger		Off Board (8A)
Output Power of Motor	hp (W)	1.1 (820)
Maximum Speed ²	mph (km/hr)	8 (13)
Maximum Climbing Angle / Safe Climbing Angle(based on user weight of 90kg/200lb)	% (°)	14 (8)
Range (per charge)(After the battery and mechanical moving parts fully break in)	mile (km)	49 (79)
Turning Radius	in (mm)	72 (1829)
Ground Clearance	in (mm)	5 (127)
Handle Bar	Type	T-Bar

SHOPRIDER

Toll-Free: 800.743.0772

www.shoprider.com