


# FLAGSHIP

TE-889XLSN

Introducing Shoprider's top of line FLAGSHIP scooter. The name says it all. The ultra exclusive FLAGSHIP will satisfy those with discriminating taste for only the very best. The FLAGSHIP's stylish canopy enclosure provides comfortable, all season transportation. Sturdy, key locking and removable doors provide security as well as convertability. The FLAGSHIP also comes standard with a fully operational lighting system, full suspension for a luxurious ride, even a windshield wiper/washer system with reservoir. Trade in the Rolls-Royce for your FLAGSHIP today. Enjoy tomorrow's products today.. from Shoprider.

## Available Colors



- **350 lb. weight capacity**
- **Front and Rear Suspension**
- **Wiper and rear view mirrors are standard equip**
- **Genuine warm weather safety windshield automotive**
- **Windshield wiper & automatic wiper fluid**

Overall Dimensions (L1 x W x H)	in	63 x 29 x 64
Number, Size of Tires (front, rear)	Front (Rear)	2, 12" (2, 12")
Suggested User Weight (on a level road)	lb (kg)	Max. 350 (160)
Weight of Heaviest Part (w/o baskets and seat)	lb (kg)	48 (106)
Total Weight (w/ batteries, basket, seat)	lb (kg)	149 (328)
Each Battery Module Weight	lb (kg)	N/A
Battery Module Capacity	_V_Ah x Pcs	12V36/50/75Ah x 2
Charger		Off Board (3,4 or 8)
Output Power of Motor	hp (W)	2 (1492)
Maximum Speed <sup>2</sup>	mph (km/hr)	7(11)
Maximum Climbing Angle / Safe Climbing Angle(based on user weight of 90kg/200lb)	% (°)	15 (9)*
Range (per charge)(After the battery and mechanical moving parts fully break in)	mile (km)	75, 50 (31)
Turning Radius	in (mm)	70 (1780)
Ground Clearance	in (mm)	N/A
Handle Bar	Type	Delta

**SHOPRIDER**

Toll-Free: 800.743.0772

[www.shoprider.com](http://www.shoprider.com)