

Jazzy[®]

SELECT 14 SERIES

THE EASIEST POWER CHAIR TO USE. #1 IN MANEUVERABILITY.

The Jazzy Select 14 Series delivers the advantages of patented mid-wheel drive for maneuvering in tight, compact spaces. Both the Jazzy Select 14 and Jazzy Select 14 XL are loaded with performance and convenience features like a powerful 60 amp PG VR2 controller, larger batteries for greater per-charge range, and on-board battery charger.



Jazzy Select 14

- 300 lbs. weight capacity
- Standard 18"W x 18"D high-back contoured seat



Jazzy Select 14 XL

- 350 lbs. weight capacity
- Contoured seating available up to 22"W x 20"D



Pride[®]
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

JZ SEL 14 & 14XL SS 01/07



The Jazzy® Select 14 series models include:

- Jazzy Select 14-C
- Jazzy Select 14-SS
- Jazzy Select 14 XL-C
- Jazzy Select 14 XL-SS

FEATURES

- Patented mid-wheel drive design for tight space maneuverability
- 14" knobby tires and optimum weight distribution provide superior traction
- NF-22 batteries for maximum range
- New 60 amp, PG VR2 controller
- On-board, 4 amp charger

OPTIONS

- Contoured high-back seat
- Contoured high-back with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Safety flag


The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %.

	Jazzy® Select 14	Jazzy® Select 14 XL
Weight capacity:	300 lbs.	350 lbs.
Maximum speed:	Up to 4.5 mph	Up to 4.5 mph
Ground clearance:	3"	3"
Turning radius:	21"	21"
Overall length:	38" without front riggings	38" without front riggings
Overall width:	24.75"	24.75"
Wheels:	Drive: 14"	Drive: 14"
Front anti-tips:	6" solid	6" solid
Rear casters:	8" solid	8" solid
Suspension:	Limited	
Drive train:	Two-motor, mid-wheel drive	
Braking system:	Intelligent braking (electronic, regenerative disc brakes)	
Std. Electronics:	PG VR2, 60 amp	
Specialty controls:	N/A	
Battery charger:	On-board, 4 amp	
Per-charge range:	Up to 15 miles	
Battery requirements:	(2) 12 volt, deep cycle	
Battery size:	NF-22	
Transport weight		
Battery weight:	38 lbs. each	38 lbs. each
Base weight:	105 lbs.	106.5 lbs.
Std. Seat dimensions:	18" x 18"	18" x 20"
Warranty:	Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-months limited warranty on drive motors	

Colors:



- Full Platform
- Power chair lift with straps
- Backpacker™



Power Chair Lift

These lifts are compatible with the Jazzy® Select 14 Series.



Backpacker™

Pride[®]
Mobility Products Corp.