



F O L D A B L E P O W E R W H E E L C H A I R

Roll

Alero

Foldable electric Power Chair available in three sizes.



Rise.

Roll.

Rest.

Ride.



YOUR LIFE IN MOTION

Alero GP52 Specifications

Model:	16"	18"	20"
Weight:			
Base (w/out footrest)	84 lbs	85 lbs	87 lbs
Overall Dimensions:			
Length	47"	47"	49"
Width	24.5"	26.5"	28.5"
Load Capacity:	250 lbs	250 lbs	250 lbs
Seat Dimensions:			
Depth	16"	16"	18"
Width	16"	18"	20"
Operating Range:	20 miles*		
Top Speed:	3.7 mph		
Turning Radius:	27.5"		
Wheel Sizes:			
Drive Wheels	12.5" x 2.25" Foam-filled		
Casters	8" x 1.75" Foam-filled		
Anti-tippers	2" (Rear Mounted)		
Controller:	PG VSI-50A		
Charger:	Off board, 5 amp		
Brakes:	Electromagnetic		
Color:	Black		
Seat Color:	Black		
Battery Type	(2) U-1 36Ah		

Power Chair Options

- Elevating Leg Rests
- Specialty Leg Support
- Left or Right Mounted Controller

Features

- Flat Free Solid Tires
- Swing Away Footrests
- Clothing Guards
- Adjustable Controller Bar (length)
- Adjustable Rear Anti-tip Wheels
- Foldable Battery Tray
- Strong Plastic Battery Containers
- Quick Plug Battery Connections
- Charger Pocket
- Seat Belt
- Folding Backrest
- Seat Cushion
- Seat Loop Handle for Easy Folding



* Battery range will vary due to rider weight, drive surface, drive slope, and tire pressure.

Literature is current at time of printing. Company reserves the right to make changes.