

P O W E R C H A I R



Seat Color



Candy Apple Red



Onyx



Body Colors

Blue



Hunter Green



**New Standard Design**

# Roll

golden  
**COMPASS**<sup>™</sup>

True, stable platform, center-wheel drive power chair with an articulating chassis for the ultimate in maneuverability and stability.



Rise.

Roll.

Rest.

Ride.



**YOUR LIFE IN MOTION**

## Golden Compass™ Specifications

<b>Weight:</b>	
Base	88 lbs
Seat	25 lbs
Total Weight	176 lbs (with U1 batteries)
<b>Overall Dimensions:</b>	
Length	38.5"
Width	24.5"
Seat Height from ground	18" (with U1 batteries) adjustable to 23"
<b>Top Speed:</b>	5 mph
<b>Operating Range:</b>	20 miles*
<b>Load Capacity:</b>	300 lbs
<b>Turning Radius:</b>	19"
<b>Wheel Sizes:</b>	
Drive Wheels	10" x 3.5"
Casters	6" x 2"
<b>Seat Dimensions:</b> Standard: 18" x 16"	
These sizes available at an upcharge	
Max. dimension	18" x 18", 18" x 20", 20" x 18"
Min. dimension	24" x 20"
	16" x 16"
<b>Footrest:</b>	
Height adjustable	Yes
Angle adjustable	Yes
Length adjustable	Yes
<b>Controller:</b> GP-600 Standard: Dynamic A-Series 50 amp, fully programmable integral	
Optional	PG VR2, 50 amp, programmable Off Board Charger
<b>Charger:</b>	
<b>Colors:</b> Candy Apple Red, Blue, Onyx, Hunter Green	
<b>Seat Colors:</b> Standard: Black Vinyl, Stadium Style	
Optional Two-tone:	Coffee/Sand, Charcoal/Gray Captain's
<b>Battery Type</b> (recommended):	
	2 (two) 12 Volt U-1 or
	2 (two) 12 Volt 22NF



Optional seating featuring our exclusive fold-flat seat hinge design for easy transport

## Power Chair Accessories and Options

- Captain's seat
- High back seat
- Pro Flex seat
- Van/pan seat
- Swing-away legrests
- Elevating legrests
- Stump support
- Footplate extender
- Positioning belt
- Cane/crutch holder kit
- Walker holder kit
- Oxygen tank holder kit
- Pack N Go

## 4 exciting colors



Candy Apple Red



Blue



Hunter Green



Onyx

\* Battery range will vary due to rider weight, drive surface, drive slope.  
Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.