

AmeriGlide

AmeriGlide 300



The AmeriGlide 300 easily fits into most vehicles and is a perfect complement to scooters and power chairs from Pride and other manufacturers. Engineered to provide smooth, efficient operation, maximum durability, and sporty style, the AmeriGlide 300 makes it easy and affordable to take your mobility product with you wherever you go.

- Lifting capacity 300 lbs.
- Low headroom required - 30"
- Telescoping base legs and extending boom arm length for maximum installation adjustability
- Boom arm end cap to protect vehicle headliners
- End swivel hook to prevent strap twisting
- Stop warning on strap to identify lifting limit
- Lift is easily disassembled for stowing
- Integrated boom arm securement when not stowed
- Remote hand control operation
- Integrated dual protective devices to prevent lift or vehicle electrical system damage
- Interchangeable appearance panel to complement your vehicle's interior
- *COMING SOON!* Optional power rotating base
- *COMING SOON!* Wireless remote operation

AmeriGlide

www.ameriglide.com

AmeriGlide 300



Can easily lift the following Jazzy® Power Chairs ¹

- Jazzy® 600
- Jazzy® 610
- Jazzy® 1103 Ultra
- Jazzy® 1107
- Jazzy® 1113 ATS
- Jazzy® 1121
- Jazzy® 1122
- Jazzy® Select
- Jazzy® Select 14
- Jazzy® Select 14 XL

Can easily lift the following Quantum® Power Chairs ¹

- Dynamo® 1107
- Quantum® 600
- Quantum® 610
- Quantum® 1107
- Quantum® 1121
- Quantum® 1402
- Quantum® 1420
- Vibe®

Can easily lift the following Pride® Scooters

- Celebrity® X
- Go-Go Elite Traveller™
- Go-Go® Ultra X
- Hurricane®
- Legend®
- Legend® XL
- Maxima
- Rally®
- Revo®
- Sonic®
- Sundancer®
- Victory®

Can easily lift the following Jet™ Power Chairs ¹

- Jet™ 3 Ultra
- Jet™ 7
- Jet™ 10 Ultra

¹ Seating system dependent

Notes: The AmeriGlide 300 not only works great with Pride Mobility products, but also with other power chair and scooter brands.