




Rally® 3-wheel

Rally® 4-wheel

Model number:	SC151	SC155
Weight capacity:	300 lbs.	300 lbs.
Maximum speed:	Up to 4.25 mph	Up to 4.25 mph
Ground clearance:	3"	3"
Turning radius:	39"	54"
Overall length:	41"	44"
Overall width:	21"	21"
Wheels:		
Front:	2" x 8" pneumatic	2" x 8" pneumatic
Rear:	2.5" x 9" pneumatic	2.5" x 9" pneumatic
Anti-tips:	Rear	Rear
Suspension:	N/A	N/A
Transport weight		
Total w/o batteries:	93 lbs.	116 lbs.
Heaviest piece when disassembled:	45 lbs. (rear section)	45 lbs. (front or rear section)
Drive train:	Sealed transaxle, 24 volt DC motor	
Dual braking system:	Electronic, regenerative and electromechanical (electronic, regenerative disc brakes)	
Power elevating seat:	N/A	
Std. Seating:	Foldable molded plastic	
Battery charger:	On-board, 3 amp (standard)	
Per-charge range:	Up to 25 miles	
Battery requirements:	(2) 12V deep cycle (AGM or gel cell type recommended)	
Battery size:	U-1 or optional 17AH	
Warranty:	3-year limited warranty	



Benefits

- Controls at your fingertips
- Easy-to-use external freewheel lever
- Aluminum alloy mags for great style
- Front basket to hold all your belongings
- Seat offers height adjustment and swivel with flip-up armrests for comfortable positioning
- Adjustable tiller for maximum comfort
- On-board accessory bracket
- Easily accessible charger port and circuit breaker
- Easily disassembles for transport

Options

- Forearm crutch holder
- Lap belt
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Heavy-duty dust cover
- Saddlebag
- Quad cane holder
- Clip, cane attachment
- Flat-free tires



- Backpacker
- Scooter double entry
- Full Platform



Double entry



Backpacker

These lifts are compatible with the Rally®.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

Pride®
Mobility Products Corp.

RALLY®

**STYLISH DESIGN & PREMIUM
PERFORMANCE AT AN
EXCELLENT VALUE.**

RALLY®

The Rally® is sure to cause a stir with its futuristic styling, technologically advanced performance and overall great value—a combination irresistible and unmatched in the small scooter class. With the Rally, you'll stake your claim on a lifestyle of fun and freedom.



Pride®
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

RALLYSS06-06