

M A X I C O M F O R T ™ S E R I E S

Shown in Night Navy Brisa®



- Mocha Suede 
- Cactus Suede 
- Burgundy Vinyl 
- Galaxy Vinyl 

New!

Rise



Milan
with
Maxi Comfort™

The latest in lift chair styling from Golden. The new Milan is the newest model in the MaxiComfort™ series. Featuring a contemporary, slim design, the Milan features a luxurious seam back and soft seat. With unique butterfly arm design, the Milan will effortlessly blend with any contemporary decor.



Rise. Roll. Rest. Ride.



YOUR LIFE IN MOTION

Milan

M A X I C O M F O R T ™ S E R I E S

Model Number	PR-607
Size	Medium
Weight Capacity	375 LB
Recline Positions	Maxi Comfort
Standard Fabrics	4
Fabric Options**	28
Chaise Seat	YES
Back Type	SEAM
Head Rest Cover	YES
Arm Covers	NO
Battery Back Up	YES
Overall Width	30"
Overall Height	43.5"
Width Between Arms	20"
Floor to Top of Seat	19.25"
Seat Depth	20.5"
Seat to Top of Back	29.5"
Seat to Top of Arm	7.5"
Floor to Back Edge of Seat Raised	26.25"
Floor to Front Edge of Seat Raised	21.75"
Seat Angle Raised	20°
Floor to Top of Foot Rest Reclined	21.5"
Floor to Top of Back Reclined	27.5"
Back Angle Closed	93°
Back Angle Reclined	162°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	180°
Distance Required from Wall Reclined	25"
Extension reclined from Back of Seat to end of Foot Rest	41"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical parts	3 YRS.
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 YRS.
Warranty on Electrical Labor	1 YR.
Removable Back	YES
Customizable *	NO
Left Hand Control Available	YES
Extra Pocket Available	YES
Heat & Massage Option Available	YES
Foot Rest extension Option Available	NO
Smart Tek System	YES

Hand-blown, Overstuffed, Conjugated Microfiber Dacron Fill
SUPER SUPPORTIVE

Northern Maple wood
• Double-doweled at Stress Points
• Reinforced Miter Blocks
• Stapled and Glued
STRONGEST IN THE INDUSTRY

2" of Block Foam
EXTRA COMFORT

Quick Release Hand Control
EASY DISCONNECT


Rear Landing Legs
EXTRA STABILITY

5 Heavy Gauge Springs
LONG LASTING FOR INCREDIBLE COMFORT

Widest Fabric Selection in Industry
FITS IN ANY HOME DÉCOR

10 Gauge Steel Scissor Mechanism
STRONG & LONG LASTING

Adjustable Anti-Skid/Scrape Feet
STABILITY & SUPPORT



501 Medium shown



501 M shown



501 M shown



501 M shown

Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.

Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.

Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.

* See Customize Chart for Dimension limits
Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.
** Value Series fabrics not available on PR607.



Smart Tek

Our luxury lift & recliners are equipped with electronics able to assure a secure lift every time you need it.



Maxi Comfort™



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Our chairs meet some of the toughest safety standards ever. Our entire chair meets ETL (Electrical Testing Lab) standards to ensure electrical safety. ETL is a nationally recognized, third party testing facility.

