

# Dasher<sup>®</sup>

(GK83, GK8)

SHOPRIDER<sup>®</sup>



Dasher<sup>®</sup> 3 (GK83)

The new **Dasher<sup>®</sup>** scooter from *Shoprider* sets another gold standard in portable mobility. Packed with industry first features, available in 3 and 4 wheel models, the **Dasher<sup>®</sup>** is the first scooter with Splash-n-Dash<sup>™</sup> quick color change moldings. It's the first completely connectorless portable scooter, with its heaviest component an industry leading 26lbs! Talk about ULTRA value...**Dashers<sup>®</sup>** have solid tires, a comfortable swivel seat with arms, and a basket included in the price without the need for an upgrade! The brand you know. The quality you trust. ULTRA portable. ULTRA lightweight. ULTRA value.



Dasher® 4 (GK8)



The **Dasher®** breaks down easily for quick transport and very user-friendly re-assembly

## Dasher® (GK83, GK8)

- Splash 'n Dash™ Quick Color Change System
- Red, Blue and Champagne Color Panels Included
- Heaviest Part Weighs Only 26 Pounds (Dasher™ 3)
- The Easiest Disassembly and Re-assembly in the Industry Virtually Eliminating the Need and Expense for Special Lifts or Ramps
- Convenient Swiveling Seat Feature While Seated
- Available in 3-Wheel and 4-Wheel Models

### Unit Specifications

		Dasher® 3	Dasher® 4
Overall Dimensions (L' x W x H)	in	42.5 x 22 x 35	44 x 22 x 33
Number, Size of Tires	Front	1, 8"	2, 7"
	Rear	2, 8"	2, 8"
Suggested User Weight (on the Level Road)	lb	Max. 250	Max. 250
Weight of Heaviest Part (Front Chassis Assembly)	lb	26	33
Total Weight (w/ 12Ah Battery Pack, Basket, Seat)	lb	89	94
Battery Pack Module Weight	lb	22	22
Battery Module Capacity	_V_Ah x Pcs	12V12Ah x 2	12V12Ah x 2
Charger		Off-Board (2A)	Off-Board (2A)
Output Power of Motor	hp	.5	.5
Maximum Speed <sup>2</sup>	mph	3.8	3.8
Safe Climbing Angle (based on user weight of 100kg/220lb) / Max. Climbing Angle (based on user weight of 75kg/165lb)	% (°)	10 (6) / 21 (12)	10 (6) / 21 (12)
Range <sup>3</sup> (per charge w/ STD 12AH battery / 90kg User Weight) (After the battery and mechanical moving parts fully break in)	mile	7.5	7.5
Turning Radius	in	33	37
Ground Clearance	in	1.5	1.5
Handle Bar	Type	T-bar	T-bar

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

- 1) Include the anti-tip wheel or the rear castor.
- 2) Driver weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown here:

- a) The weight of occupant
- b) Ground surface
- c) Battery capacity and conditions
- d) Type of charger
- e) Ambient temperature

- f) The way of driving
- g) If the battery and mechanical moving parts fully break in
- h) Etc.
- i) Height is measured at the lowest point



Connector-less Technology



Convenient Built-in Charger Compartment



Adjustable Tiller

# SHOPRIDER®

Toll-Free: 800.743.0772

www.shoprider.com