

SHOPRIDER®

Shopperider Mobility Products, Inc.



GOES EVERYWHERE YOU GO

LIGHTWEIGHT-Only 112 lbs total weight

PORTABLE-Small & manageable – fits in the car trunk

- ◆ Durable & Portable
- ◆ Forward & Reverse Controls
- ◆ Variable Speed Settings
- ◆ Adjustable Seat, Handlebar & Footrest for Maximum Comfort
- ◆ Battery-Life Indicator
- ◆ Automatic Electromagnetic Braking System
- ◆ Anti-Tip Wheels for Added Safety
- ◆ Horn & Reverse Alarm
- ◆ Quick-Release Battery-Pack
- ◆ Charge Battery-Pack without Unit
- ◆ Flat Free Polyurethane Tires

World Patented
Cutie "L" SCOOTER

SPECIFICATIONS	
Overall Dimensions LxWxH	43" x 26" x 36"
Size of Tires	Front 9" x 1, Rear 11" x 2
Total Unit Weight	112 lbs
User Weight Capacity	300 lbs
Range (per charge)	20 miles
Battery Capacity	12v 21Ah x 2
Charger	3 amp OFF-BOARD
Turning Radius	35"
Ground Clearance	3"
Maximum Speed	5 mph
Drive Train	Sealed Transaxle, 24v DC motor