

The Tush Push

Tush Push

The Tush Push is the most exciting aid for lifting and lowering a person to a toilet. If you have lost the body strength to get up from the toilet, then our Tush Push is the right solution for you.



Our customers say their lives have greatly improved after buying the Tush Push.



Seat lifts for male user.

Features

- ✓ Safe, Strong, Dependable and Portable.
- ✓ Easy to install on most toilets.
- ✓ Adjustable and removable arms.
- ✓ Touch control handset for lifting action.
- ✓ Lifts to natural body movement.



Remember- If you have a lift chair in the living room, you need a Tush Push in the bathroom.



- ✓ Removable plastic shield for cleaning.

3 Year Limited Warranty

Technical



Power Supply:	115vac to 24vdc
Weight:	34 lbs
Weight Capacity:	275 lbs.
Seat Width:	14 in.
Total Width:	24 in.
Seat Height min:	19 in.
Seat Height max:	26 in.
Arm Height max:	36 in.



2 Motor option for 500lbs lifting power.