



Featuring new 2 button lighted hand control.

NEW!

ELEGANCE COLLECTION LL-510 SERIES

Pride® Lift Chairs | *Experience the Luxury*

The LL-510 3-Position Chaise Lounger from Pride delivers unique style, smooth operation and maximum durability. The LL-510 features a contemporary contoured pillow back for enhanced lumbar support and exceptional comfort.

BENEFITS

- New look of the Elegance Collection series
- Contemporary contour pillow back design
- Available in 3 sizes S, M & L
- Over 60 fabric choices to select from
- Additional lumbar support
- 375 lbs. weight capacity

	LL-510S	LL-510M	LL-510L
Seat To Floor:	19"	20"	21"
Seat Depth:	19"	20"	22"
Seat Width:	19"	20"	23"
Top of Back to Seat:	25"	27"	29"
Overall Chair Width:	34"	35"	37"



1-800-800-8586 (U.S.) • 1-888-570-1113 (Canada)
www.pridemobility.com • www.prideliftchairs.com

