

# S-Ergo 106

## Lightweight Ergonomic Manual Wheelchair W/ Angle Adjustable Back



- Medicare code K0004
- High strength, weighs only 27 lbs. (w/ footrests)
- AEGIS Anti-Bacterial Upholstery, removable & washable
- Has 3 easy angle adjustments: “Working” Angle (95°), “Moving” Angle (100°) & “Resting” Angle (105°)
- “Rigidifying” push-bar & folds to strengthen back support, comes standard on this model!
- Folding precision backrest
- Fixed wider/concaved armrest & Ergonomic footrest
- “Tube-in Center” foot-plate, assures side leg support
- 24” Rear polyurethane/high profile/flat free wheels
- 18x17 Seat width only
- Anti-Tippers standard on this model!
- Frame Color: Pearl Silver Only
- Removable Upper Upholstery: Silver Anti-Biotic
- Bottom Upholstery: Black breathable strong mesh
- Weight Capacity of 250 lbs

### Features and Uses For S-Ergonomic Series

- Ergonomic, lightweight, compact
- Pressure relief, stabilization, and anti-slipage
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider arm pad, providing comfort/support



Working angle(95°) helps relieving exhaustion and replenishing energy.

Moving angle(100°) allows effort-saving flexibility and dynamic display.

Resting angle(105°) offers relaxation of body and mind, and recharging of vitality.

items	F/R wheels												Max. User Weight
-------	------------	--	--	--	--	--	--	--	--	--	--	--	------------------

S-Ergo106F16 ~ 16"x17"	7"/24"	18"	17"	23"	12"	9"	39"	19"	17 ½"	35 ½"	27 lbs.	250 lbs.
S-Ergo106F18 ~ 18"x17"	7"/24"	18"	17"	25"	12"	9"	39"	19"	17 ½"	35 ½"	27 lbs.	250 lbs.

**\*Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery**  
**\*\*All models come equipped with 7"x1" front polyurethane wheels**



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue. City of Industry. CA 91748-1418

Tel.626-581-2235 Fax.626-581-2335