



KM-AT60

Ultra Lightweight Active Indoor/Outdoor Wheelchair

Features:

- High-strength rigid aircraft grade alloy frame with foldable backrest providing easy carry & transportation
- Retightening straps to avoid spinal problems with seat & backrest. Also doubles for bucket adjustments!
- Folding Backrest for transporting
- Standard: Spoke type rear wheel with Quick Release function (24") & alloy push rim. Also Front 5" Polyurethane Caster
- Adjustable Back Angle
- Flip-back side panels
- Double Adjustable Handbrakes
- Only 23 lbs with a weight capacity of 220 lbs

The compact design and elegant style of the AT-60 bring you outstanding activity with allure and style!



Adjustable Seat & Back Straps



Folding Backrest



Flip-Back Side Panels

****Available Options ~ Anti-Tipper, Seat Belt, Leg-tube Protector, Magnesium Alloy Wheel & Handbrakes***



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue. City of Industry. CA 91748-1418

Tel.626-581-2235 Fax.626-581-2335