



S-Ergo 115-TP

Standard Lightweight Ergonomic Manual Wheelchair

- High strength, weighs only 22 lbs. (w/o footrests)
- Anti-Bacterial Upholstery
- Folding backrest
- Fixed armrest, Swing away footrest
- “Tube-in Center” foot-plate, assures better side leg support
- Extending brake handles
- 14” Rear polyurethane, high profile, flat free wheels
- 16x17 or 18x17 Seat width
- Frame Color: Pearl Silver & Rose Red
- Upholstery: Silver/Black
- Weight Capacity of 250 lbs

Features and Uses for S-Ergo Series

- Ergonomic, lightweight, compact
- Pressure relief, stabilization, and anti-slipage
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider arm pad, providing comfort/support

Items	F/R wheels											Max. User Weight
Width												

S-Ergo 115 ~ 16”x17”	7”/24”	16”	17”	23”	12”	8”	35 ¾”	18 ½”	16”	36”	25 lbs.	250 lbs.
S-Ergo 115 ~ 18”x17”	7”/24”	18”	17”	25”	12”	8”	35 ¾”	18 ½”	16”	36”	25 lbs.	250 lbs.

**Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery, Canopy, Rigidifying Push-Bar*

***All models come equipped with 7”x1” front polyurethane wheels*



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue. City of Industry. CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335