

SPECIALTY COLLECTION



Small
LC-770S^{4,12}
Infinite-Position, Trendelenburg
Chaise Lounger
5'4" & Below



Large
LC-770L^{3,4,12}
Infinite-Position, Trendelenburg
Chaise Lounger
5'5" - 6'0"



Medium
LC-570T¹²
3-Position
Chaise Lounger w/ Transfer Arms
5'5" - 6'0"



Extra-tall
LC-575¹²
3-Position, Chaise Lounger
w/ Footrest Extension
6'1" - 6'6"



Large
LC-571¹²
3-Position
Chaise Lounger w/ Flip-open Arms
5'5" - 6'0"



Medium
LC-805^{4,7,12}
2-Position
Chaise Lounger Wall Hugger
5'4" - 6'0"



Petite Wide
LC-525PW⁸
Infinite-Position "Zero-gravity"
Trendelenburg Chaise Lounger
5'2" & Below



Small
LC-525S⁸
Infinite-Position "Zero-gravity"
Trendelenburg Chaise Lounger
5'3" & Below



Medium
LC-525M⁸
Infinite-Position "Zero-gravity"
Trendelenburg Chaise Lounger
5'4" - 6'0"



Large
LC-525L⁸
Infinite-Position "Zero-gravity"
Trendelenburg Chaise Lounger
5'9" - 6'2"



Super Petite
LC-102⁴
3-Position Chaise Lounger
4'8" - 5'4"

AVAILABLE ON ALL LIFT CHAIRS, EXCEPT THE LC-102 & LC-105

NEW CLOUD 9 STANDARD FABRICS




Medium
LC-105⁴
3-Position Chaise Lounger
5'2" - 5'10"

LIFT CHAIR SPECIFICATIONS

AMERICA'S #1
LIFT CHAIR BRAND



	Weight Capacity	Seat-to-Floor	Seat Depth	Seat Width	Top of Back to Seat	Back Style	Fully Padded Chaise	Heat & Massage Option	Total Comfort Seating*	Distance From Wall	Overall Chair Width	Chair Weight	
Elegance Collection	LC-570S	375 lbs.	19"	19"	19"	24.5"	Blown Fiber Pillow Back	Yes	Yes	Yes	16.5"	35"	110 lbs.
	LC-570M	375 lbs.	20"	20"	20"	26.5"	Blown Fiber Pillow Back	Yes	Yes	Yes	19"	36"	113 lbs.
	LC-570W	375 lbs.	20"	20"	24"	26.5"	Blown Fiber Pillow Back	Yes	Yes	Yes	19"	40"	117 lbs.
	LC-570L	375 lbs.	21"	22"	22"	27"	Blown Fiber Pillow Back	Yes	Yes	Yes	24"	38"	114 lbs.
	LC-550M	375 lbs.	19"	20"	21"	27"	Blown Fiber Split Back	Yes	Yes	No	19"	38"	119 lbs.
	LC-585	375 lbs.	20"	20"	20"	26"	Blown Fiber Biscuit Back	Yes	Yes	No	18"	36"	120 lbs.
	LC-450	375 lbs.	20"	20"	20"	27"	Blown Fiber Split Back	Yes	Yes	No	17"	38"	119 lbs.
Specifications on all Elegance Collection Lift Chairs:		Wood Frame: Engineered Furniture Grade Laminate/Hardwood Frame, Warranty: 7-Year Prorated, Lifetime Limited on lift mechanism, shell and recline mechanism.											
Specialty Collection	LC-770S	375 lbs.	19"	19.5"	19"	24"	Sewn Pillow Back	Yes	Yes	No	28"	34"	129 lbs.
	LC-770L	375 lbs.	21"	21"	22"	28"	Sewn Pillow Back	Yes	Yes	No	30"	37"	136 lbs.
	LC-570T	375 lbs.	20"	20"	20"	29.5"	Blown Fiber Pillow Back	Yes	Yes	No	20.5"	35.5"	128 lbs.
	LC-575	375 lbs.	21"	22"	22"	31"	Blown Fiber Pillow Back	Yes	Yes	No	24"	41.5"	138 lbs.
	LC-571	375 lbs.	21"	22"	22"	27"	Blown Fiber Pillow Back	Yes	Yes	No	16"	38"	133 lbs.
	LC-805	300 lbs.	19"	21"	20"	27"	Blown Fiber Sewn Split Back	Yes	Yes	No	4"	35"	118 lbs.
	LC-525PW	375 lbs.	19"	19"	24"	25.5"	Blown Fiber Split Back	Yes	No	No	28"	38"	130 lbs.
	LC-525S	375 lbs.	19"	18"	18"	25.5"	Blown Fiber Split Back	Yes	No	No	28"	32"	123 lbs.
	LC-525M	375 lbs.	20"	20"	20"	26.5"	Blown Fiber Split Back	Yes	No	No	29"	33"	126 lbs.
	LC-525L	375 lbs.	21"	22"	22"	27.5"	Blown Fiber Split Back	Yes	No	No	30"	35"	128 lbs.
	LC-102	250 lbs.	16"	18"	16"	23.5"	Dual Pillow Back	Yes	No	No	16"	25.5"	96 lbs.
	LC-105	325 lbs.	18"	19"	20"	27.5"	Sewn Pillow Back	Yes	No	No	18"	32.75"	96 lbs.
Specifications on all Specialty Collection Lift Chairs:		Wood Frame: Engineered Furniture Grade Laminate/Hardwood Frame, Warranty: 7-Year Prorated, Lifetime Limited on lift mechanism, shell and recline mechanism.											
Heritage Collection	LC-358P	375 lbs.	18"	18"	18"	24"	Blown Fiber Button Back	Yes	Yes	Yes	17.5"	34.5"	112 lbs.
	LC-358PW	375 lbs.	19"	19"	24"	23.5"	Blown Fiber Button Back	Yes	Yes	Yes	17"	39"	118 lbs.
	LC-358S	375 lbs.	19"	19"	18"	24"	Blown Fiber Button Back	Yes	Yes	Yes	17"	35"	118 lbs.
	LC-358M	375 lbs.	20"	20"	20"	24.5"	Blown Fiber Button Back	Yes	Yes	Yes	17"	35"	117 lbs.
	LC-358L	375 lbs.	21"	22"	22"	25.5"	Blown Fiber Button Back	Yes	Yes	Yes	17.5"	37"	120 lbs.
	LC-358XL	500 lbs.	20"	22"	26"	29"	Blown Fiber Button Back	Yes	Yes	Yes	18"	39.5"	159 lbs.
	LC-358XXL	600 lbs.	20"	22"	30"	29"	Blown Fiber Button Back	Yes	Yes	Yes	18"	43"	168 lbs.
Specifications on all Heritage Collection Lift Chairs:		Wood Frame: Engineered Furniture Grade Laminate/Hardwood Frame, Warranty: 7-year Prorated, Lifetime Limited on lift mechanism, shell and recline mechanism.											
Classic Collection	LC-10	325 lbs.	20"	20"	20"	25"	Blown Fiber Button Back	No	No	No	12"	33.5"	110 lbs.
	LC-15	325 lbs.	20"	20"	20"	27"	Blown Fiber Button Back	No	No	No	18"	33.5"	119 lbs.
	LC-20	375 lbs.	20"	20"	20"	24.5"	Blown Fiber Sewn Pillow Back	No	Yes	No	11"	32.5"	114 lbs.
	LC-30	375 lbs.	20"	20"	20"	26"	Blown Fiber Split Back	Yes	Yes	No	18"	34"	118 lbs.
Specifications on all Classic Collection Lift Chairs:		Wood Frame: Engineered Furniture Grade Laminate/Hardwood Frame, Warranty: 2-year Prorated, Lifetime Limited on lift mechanism, shell and recline mechanism.											



Protect your investment with the Pride® Lift Chair Care Kit




AMERICA'S #1 LIFT CHAIR BRAND

Pride® Lift Chairs

LIFT CHAIR GUIDE

www.pridemobility.com • www.prideliftchairs.com

OTHER AVAILABLE FABRICS



ELEGANCE COLLECTION



Medium
LC-550M^{8,12}
3-Position
Chaise Lounger
5'4" - 6'0"



Medium
LC-585^{6,12}
3-Position
Chaise Lounger
5'5" - 6'0"



Medium
LC-450^{6,12}
3-Position Chaise Lounger
5'4" - 6'0"



Small
LC-570S^{8,12}
3-Position
Chaise Lounger
5'3" & Below



Medium
LC-570M^{8,12}
3-Position
Chaise Lounger
5'4" - 6'0"



Medium Wide
LC-570W^{8,12}
3-Position
Chaise Lounger
5'4" - 6'0"



Large
LC-570L^{8,12}
3-Position
Chaise Lounger
5'9" - 6'2"



Medium
LC-358M^{6,8,12}
3-Position
Chaise Lounger
5'4" - 6'0"



Large
LC-358L^{6,8,12}
3-Position
Chaise Lounger
5'9" - 6'2"



Petite Wide
LC-358PW^{8,12}
3-Position
Chaise Lounger
5'3" & Below



Small
LC-358S^{8,12}
3-Position
Chaise Lounger
5'3" & Below



Heavy-duty
LC-358XL^{1,4,7,10,12}
3-Position
Chaise Lounger
5'6" - 6'1"



Very Heavy-duty
LC-358XXL^{1,4,7,10,12}
3-Position
Chaise Lounger
5'6" - 6'1"

HERITAGE COLLECTION

AVAILABLE ON ALL LIFT CHAIRS, EXCEPT THE LC-102 & LC-105

NEW CLOUD 9 STANDARD FABRICS



BLACK CHERRY



PACIFIC



SPRUCE



STONE



WALNUT

Classic COLLECTION



Medium
LC-10⁸
2-Position
5'2" - 5'10"



Medium
LC-15⁶
3-Position
5'2" - 5'10"



Medium
LC-20
2-Position
5'4" - 6'0"



Medium
LC-30^{6,11}
3-Position
Chaise Lounger
5'4" - 6'0"

AVAILABLE ON ALL LIFT CHAIRS, EXCEPT THE LC-102 & LC-105

NEW CLOUD 9 STANDARD FABRICS



BLACK CHERRY



PACIFIC



SPRUCE



STONE



WALNUT

AMERICA'S #1
LIFT CHAIR BRAND



Look for the Quick Ship Indicator for the shortest lead times available.

Most Elegance, Specialty, Heritage and Classic Collection Lift Chairs come standard with head and arm covers. Exceptions include the LC-310 and LC-770 Series, which only come standard with a head cover; and the LC-105, which do not come standard with head or arm covers.

All lift chairs have integrated emergency battery backup.² Transformer features self-diagnostic electronics with 12' cord.

A 2-position is a TV recline.

A 3-position is a full recline.

An infinite-position is a sleep recline.

Trendelenburg position elevates feet higher than the head.

"Zero-gravity" positions thighs and lower legs higher than torso.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

1. Removable, fully padded chaise.
2. Batteries not included.
3. Available with bed option (mattress dimensions: 20"W x 73"L x 3"D).
4. Quick Ship with white toggle control.
5. Crypton Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.
6. Available with footrest extension.
7. One side pocket.
8. Fiber fill can be added or removed for custom comfort adjustment.
9. Unique spring design and a one-piece soft foam seat for outstanding relaxation.
10. Features Dual Motor Sync System. Not available with heat & massage.
11. Available in CAL133 fabric.
12. Available with Deluxe Heat & Massage.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Note: All lift chairs are shipped in an impact-resistant, reinforced, corrugated box.

