

PERSONAL TRANSPORT CHAIR



Seat Colors
Black Vinyl

Body Color
Red

Roll

New!

GOLDEN
LiteRider[®]ptc

The LiteRider brand now includes power chairs! We've taken our most popular scooter design and created the most unique power chair available! The LiteRider Personal Transport Chair is a 275 lb. weight capacity, super-compact, rear-wheel drive power chair for the consumer looking for high quality at a great value. Its unique design incorporates a basket under the seat for valuables, a comfortable, full size stadium seat and non-marking black tires on steel rims with shiny hubcaps.



Rise. Roll. Rest. Ride.

GOLDEN
technologies

YOUR LIFE IN MOTION

Model Number	GP160
Weight Capacity	275 lbs.
Drive Wheels	Rear
Maximum Speed	3.5 mph
Operating Range ¹	9.5 miles
Ground Clearance	2.5"
Turning Radius	26"
Type Batteries	2-18 AH
Colors	Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	PG VR-2 50 Amp Remote
Dimensions:	
Length	33"
Width	22.5"
Height (ground to top of back of seat)	33" - 34"
Height (ground to top of back of headrest)	N/A
Ground to Top of Seat	20.5" - 21.5"
Number of Seat Height Positions	2
Size of Increments	1.0"
Weight of Unit: (assembled)	
Base	67 lbs.
Seat w/Arms	22 lbs.
Batteries	25 lbs/pack
Tires:	
Drive Wheels	9"
Casters	6"
Anti-Tip Wheels	3"
Standard Seat:	
Slide Seat	N/A
Back Height - no headrest	16"
Back Height with headrest	N/A
Width x Depth	17" x 16"
Color	Black Vinyl
Optional Seat:	
Back Height - no headrest	16"
Back Height with headrest	N/A
Width x Depth	18" x 16", 20" x 18"
Color	Black Vinyl
Footrest:	
Height adjustable	Yes
Angle adjustable	Yes
Warranty:	
Frame	Five Years
Drive Train	13 Month Warranty
Electronics ²	13 Month Warranty

¹ Battery range will vary due to rider weight, drive surface and terrain.
² Electronics warranty excludes batteries.

Disassembles easily for transport!



1. Remove seat



2. Pull out battery pack



3. Lightweight base!



Perfect for transport!

Available Accessories:

- Cane Holder Tube
- Oxygen Tank Holder
- Cup Holder
- Quad Cane Holder
- Rear Basket
- Walker Holder
- Crutch Holder
- Safety Flag
- Power Chair Cover
- Touch-up Paint
- Deluxe Pack N' Go
- Seat Belt (Optional)

Some accessories require a mounting bracket or mounting clips.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.