



- Separates into two pieces for easy carrying
- Rotation of pin either locks or allows ramp to be separated into two pieces (shown below)
- Folds in half, carries like a suitcase
- Patented hinge eliminates pinch point
- Closure strap locks ramp panels together
- Lightweight, easy-to-handle & set up
- Durable welded construction
- Safety guide indicates if ramp is on a safe slope
- Full platform provides excellent stability
- Anti-slip, high traction surface
- Accommodates wheelchairs and scooters with various wheel configurations
- Safety DVD and steel security pins included

Listed by Underwriters Laboratories Inc.



Custom widths are available. Give your specifications to your local PVI retailer!

Item Number	Length	Width	Folded Size	Weight	Capacity*
WCR530	5'	30"	32 ¹ / ₄ " x 16" x 9"	30 lbs.	600 lbs.
WCR630	6'	30"	38 ¹ / ₄ " x 16" x 9"	35 lbs.	600 lbs.
WCR730	7'	30"	44 ¹ / ₄ " x 16" x 9"	39 lbs.	600 lbs.
WCR830	8'	30"	50 ³ / ₄ " x 16" x 9 ¹ / ₂ "	52 lbs.	600 lbs.
WC1030	10'	30"	62 ³ / ₄ " x 16" x 9 ¹ / ₂ "	63 lbs.	500 lbs.
WC1230	12'	30"	74 ³ / ₄ " x 16" x 9 ¹ / ₂ "	73 lbs.	500 lbs.

* Maximum ramp capacity: 300 lbs. (one axle) – 600 lbs. (two axles). Except on WC1030 & WC1230. Maximum ramp capacity: 250 lbs. (one axle) – 500 lbs (two axles).

