

LC-576 SERIES



The new LC-576 infinite-position, “zero-gravity,” Trendelenburg chaise lounger is the latest advancement from Pride® Lift Chairs. Featuring a waterfall back with removable filling in the back pillows for custom lumbar support, the LC-576 provides comfort and style.

Features and Benefits:

- Infinite-position
- One-button hand control push to “zero gravity” position
- True Trendelenburg position
- Removable filling in back pillows for custom lumbar support
- Industry exclusive dual pockets for convenience
- Available in Medium and Large

LC-576M

Infinite-Position
“Zero Gravity”
Trendelenburg
Chaise Lounger
5'4" - 6'0"

LC-576L

Infinite-Position
“Zero Gravity”
Trendelenburg
Chaise Lounger
5'9" - 6'2"

Weight Capacity	375 lbs.	375 lbs.
Seat-to-Floor	20"	21"
Seat Depth	20"	22"
Seat Width	20"	22"
Top of Back to Seat	26.5"	27.5"
Back Style	Blown Fiber Waterfall Back	Blown Fiber Waterfall Back
Fully Padded Chaise	Yes	Yes
Distance From Wall	29"	30"
Overall Chair Width	33"	35"
Chair Weight	126 lbs.	128 lbs.
Warranty	7-Year Prorated	7-Year Prorated

AVAILABLE FABRICS

CLOUD 9 FABRICS



BLACK
CHERRY



PACIFIC



SPRUCE



STONE



WALNUT

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.