

H E A V Y D U T Y S C O O T E R



Seat Color
Charcoal/Vinyl

Body Colors
Red
Silver

Ride

golden **Avenger**

4 wheel

The Golden Avenger four-wheel, heavy duty scooter provides an impressive blend of rugged outdoor maneuverability while giving a smooth, comfortable ride over a variety of nature's most challenging terrains. Featuring a 500 lb. weight capacity, the Avenger comes standard with NF22 batteries providing a range of 18 miles.



Rise. Roll. Rest. Ride.



YOUR LIFE IN MOTION

GA-541

Avenger 4 wheel Specifications

Model	GA541
Medicare Code	Not Classified
Weight Capacity	500 lbs
Drive Wheels	Rear
Maximum Speed	7 mph
Operating Range ¹	18 miles
Ground Clearance to Center Deck	4"
Ground Clearance to Transaxle/Motor	4"
Turning Radius	67"
Type Batteries	2-NF22
Colors	Silver, Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	PG 110 Amp
Dimensions:	
Length	57.5"
Width	24.5"
Height (ground to top of back of seat)	42.5" - 45.5"
Height (ground to top of back of headrest)	50" - 56"
Ground to Top of Seat	25" - 28"
Ground to Top of Deck	6.5"
Deck to Top of Seat	18.75" - 21.75"
Number of Seat Height Positions	4
Size of Increments	1"
Front Axle to Rear Axle	38.5"
Weight of Unit: (assembled)	268 lbs
Front Half	68 lbs
Rear Half	68 lbs
Base	N/A
Seat w/Arms	55 lbs
Batteries	37 lbs ea
Tires:	
Front	Pneumatic
Rear	13"
Casters	13"
Rear Anti-Tippers	N/A
Rear Anti-Tippers	3"
Standard Seat:	
Slide Seat	Captain's Seat
Slide Seat	N/A
Back Height - no headrest	21" - 24"
Back Height with headrest	27.5" - 31.5"
Width x Depth	20" X 17"
Color	Gray Vinyl
Warranty:	
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years

¹ Battery range will vary due to rider weight, drive surface and terrain. Foam-filled tires are optional. The operational mileage will be less if scooter is equipped with foam filled tires.


² Electronics warranty excludes batteries.



Superb Style

Features

- 500 lb weight capacity
- Speeds up to 7 mph
- Bright safety lights
- 22 NF batteries
- Up to 18 Mile Range
- 4-Way Flashers
- Great Warranty
- Oversized tires
- Oversized seat
- Removable basket



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.