

S I G N A T U R E S E R I E S

# Rise



Regal  
PR-751TY



Winston  
PR-410



Royal  
PR-752



Space Saver  
PR-906

## Signature Series

The Signature Series offers four models, each with its own unique, luxurious style. Choose from the Regal with dual open arm storage; the Winston with overstuffed, generous proportions; the Royal with its plush, three-pillow back; or the Space Saver that requires only 6" from the wall for a full recline. The Signature Series offers a style like no other.



Rise. Roll. Rest. Ride.



YOUR LIFE IN MOTION

S I G N A T U R E S E R I E S



	Winston	Regal	Royal	Space Saver
<b>Model Number</b>	<b>PR-410</b>	<b>PR-751TY</b>	<b>PR-752</b>	<b>PR-906</b>
<b>Size</b>	Large	Medium	Medium	Medium
<b>Weight Capacity</b>	375 Lb.	375 Lb.	375 Lb.	250 Lb.
<b>Recline Positions</b>	3	3	3	3
<b>Standard Fabrics</b>	4	4	4	4
<b>Fabric Options</b>	4	32	32	32
<b>Chaise Seat</b>	YES	YES	YES	YES
<b>Back Type</b>	Biscuit	3 Pillow	3 Pillow	Seam
<b>Head Rest Cover</b>	YES	YES	YES	YES
<b>Arm Covers</b>	NO	NO	YES	NO
<b>Battery Back Up</b>	YES	YES	YES	YES
<b>Overall Width</b>	43"	36"	32.5"	35.5"
<b>Overall Height</b>	46"	44"	44"	43"
<b>Width Between Arms</b>	21.5"	22"	21"	20"
<b>Floor to Top of Seat</b>	21"	20"	19.5"	22"
<b>Seat Depth</b>	21"	21.5"	21"	22"
<b>Seat to Top of Back</b>	29"	28"	28"	27"
<b>Seat to Top of Arm</b>	6.5"	8"	6"	8"
<b>Floor to Back Edge of Seat Raised</b>	29"	28"	27.5"	30"
<b>Floor to Front Edge of Seat Raised</b>	24.5"	24"	23.5"	24"
<b>Seat Angle Raised</b>	20°	18°	15°	24°
<b>Floor to Top of Foot Rest Reclined</b>	23.5"	24"	23.5"	21.5"
<b>Floor to Top of Back Reclined</b>	36"	33"	33"	37"
<b>Back Angle Closed</b>	110°	110°	115°	120°
<b>Back Angle Reclined</b>	158°	140°	145°	140°
<b>Foot Rest Angle Closed</b>	90°	90°	90°	90°
<b>Foot Rest Angle Opened</b>	180°	180°	180°	180°
<b>Distance Required from Wall Reclined</b>	19"	18"	19"	6"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	45"	42"	42"	44"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	YES	YES	YES	1 Yr.
<b>Warranty on Electrical Parts</b>	3 Yrs.	3 Yrs.	3 Yrs.	1 Yr.
<b>Pro Rated Years 4 Thru 7</b>	YES	YES	YES	NO
<b>Warranty on Mechanical Labor</b>	3 Yrs.	3 Yrs.	3 Yrs.	1 Yr.
<b>Warranty on Electrical Labor</b>	1 Yr.	1 Yr.	1 Yr.	1 Yr.
<b>Removable Back Standard</b>	YES	YES	YES	YES
<b>Smart Tek System Standard</b>	YES	YES	YES	YES
<b>Customizable</b>	NO	NO	NO	NO
<b>Left Hand Control Option Available</b>	YES	NO	YES	YES
<b>Extra Pocket Option Available</b>	YES	YES	YES	YES
<b>Heat and Massage Option Available</b>	YES	YES	YES	YES
<b>Foot Rest Extension Option Available</b>	NO	NO	YES	NO
<b>Straight Lift Option Available</b>	YES	YES	YES	NO



501 M shown

**Golden Craftsmanship**


Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

**Added Comfort and Safety**

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

**Smart Engineering**

Smooth, quiet lifting system offers a secure lift for the best performance.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

**Safety System**

Batteries not included. Safety System is designed to bring chair from recline to seated position in the event of a power outage.



**Smart Tek**

Our luxurious lift and recline chairs are equipped with electronics to help assure a secure lift every time you need it.



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68



Listing # E112417



**GOLDEN** technologies

YOUR LIFE IN MOTION

Golden Technologies • 401 Bridge Street, Old Forge PA 18518 • Tel: (800) 624-6374 • Fax: (800) 628-5165 • www.goldentech.com