

# EZ-ACCESS<sup>®</sup> a division of Homecare Products, Inc.

## Product Information Sheet

<b>Product Name:</b>	EZ-ACCESS <sup>®</sup> TOP LIP EXTENSION
<b>Stock Number:</b>	TLE
<b>Description:</b> Increase the usefulness of the TRIFOLD <sup>®</sup> or SUITCASE <sup>®</sup> Advantage Series Ramp with the Top Lip Extension (TLE). This optional extension increases the standard 3" lip to become 9", allowing the ramp to clear the rear bumper of SUVs and vans by simply attaching to the top lip of the ramp. Shipped as a set. Five-year limited warranty.	

<b>Usage Category:</b>	Advantage Series TRIFOLD <sup>®</sup> and Advantage Series SUITCASE <sup>®</sup> ramps
<b>Material:</b>	Aluminum

STOCK #	EXTENSION LENGTH	DIMENSIONS	WEIGHT
TLE	9"	8½" x 13⅞" x 2⅜"	5 lbs/SET