



KARMAN HEALTHCARE Inc.

Tel.626-581-2235 Fax.626-581-2335. 19255 San Jose Ave. City of Industry. CA 91748

Don't have room for big wheels and bulky, heavy wheelchairs? Can't walk too far? Get tired easily? Wish you had something handy and lightweight that you can use daily when you take your love one out? Introducing the perfect companion to wheelchairs: TV10B - more freedom for physically challenged people!

TV10B

Extraordinary Design!

Travel Chair / Super lightweight, compact-folding adult wheelchairs

Medicare Code E1038

- All 6061 T-6 aircraft-grade aluminum frame weighs **Only 18 pounds** including footrest!
- Quick n' easily folds to the size of a few large pizza boxes stacked together!
- Super lightweight, strong, sturdy construction, weight capacity 225 lbs
- Sleek styling, non-medical look
- Baked on enamel frame keeps maintenance to a minimum
- Folding footrest, flip-back padded armrests
- Swivel wheels are great for maneuvering in and out of tight places.
- Ultra-lightweight, micro cellular polyurethane tires for easy maneuverability
- Brake locks on rear wheels
- Removable, washable nylon fabric
- Colors available → Black (for 16" & 18" seat size) + Champagne (18" only)

Fold & Go!



Super Lightweight & Compact Folding!

The Karman TV10B travel chair is **simply the heavy weight champion of all super lightweight folding travel chairs on the market today**. Engineered to last, and affordable! It has a reinforced patented all 6061 T-6 aircraft-grade aluminum frame; the Karman travel chair offers maximum convenience, comfort, and security. Simply fold back the footrest, squeeze the quick-release lever on the push handle and fold down - it's that easy!

Fits anywhere, easy to fold and stash in your car trunk, check-in with your other luggage, or as a carry-on on airplanes or tour buses. Simple yet elegant styling doesn't look like a medical wheelchair at all. This could mean the difference between "You'd better go without me, dear" and "Let's go!"

Specification:

Seat Width (in.)	Seat Depth (in.)	Seat Height (in.)	Back Height (in.)	Overall Height (in.)	Overall Width (in.)	Overall Length (in.)	Weight Capacity (lbs.)
16	18	21	19	36 ¾	23	36	225
<u>FOLDED DIMENSIONS:</u>				Length- 28"	Width- 22 ¾"	Height- 14"	
18	18	21	19	36 ¾	23	36	225
<u>FOLDED DIMENSIONS:</u>				Length- 28"	Width- 24 ¾"	Height- 14"	
Back to Seat Angles: 95 °				Seat Angle (From Horizontal): 5 °			

