


**Revolutionary Patented New Design!**  
**Intelligent S-Shape Seating System**  
**S-105, S-115, S-305**

This intelligent S-shape Seating System provides:

- Pressure relief, stabilization, anti-slippage
- Removable upholstery is equipped with AEGIS Microbe Shield antibiotic technology, odor/staining/deterioration control
- Ergonomic, Lightweight, Compact
- Wider Arm Pad Provide Comfort Support/ Folding Backrests
- International Patents

**S-Shape Seat Frame**



**S-105**

**All 6061 T-6 Aircraft-grade Aluminum**  
**All Medicare Approved!**

- High strength, weighs **only 27 lbs** (w/foot)
- Folding backrest
- Fixed armrests & footrests
- “tube-in center” footplate assures better side support
- 16”x17 or 18”x17” seat width
- 24” rear spoke flat free wheels
- Frame Color: Pearl Silver/Rose Red, Upholstery: Silver + Black
- Weight cap 250 lbs. For all 3 models
- **Additional Options:** Companion Hand brakes




**S-115**

- High strength, weighs **only 25 lbs.** (w/o footrest)
- Folding backrest
- Removable & washable Upholstery
- 16”x17, 18”x17” or 20”x17” seat width
- Fixed Armrests, Swing-away footrests
- 24” rear spoke flat free wheels
- Frame Color: Pearl Silver/Rose Red, Upholstery: Silver + Black
- **Additional Options:** Quick Release Rear Wheels & 20” Wide Seat



**S-305**

- Light weighs **only 29 lbs.** (w/o ft)
- Folding backrest
- 16”x17 or 18”x17” seat width
- Height adjustable armrests
- Easy release lever for footrests
- Inward/outward swing-away footrests for easy operation
- Seat height adjustable
- Quick release rear spoke flat free wheels
- Frame Color: Pearl Silver/Rose Red, Upholstery: Silver + Black
- 



Karman Healthcare Inc.

19255 San Jose Ave., City of Industry, CA 91748

Tel:(626) 581-2235 Fax: (626) 581-2335

info@karmanhealthcare.com---- www.karmanhealthcare.com