



#### *Features and Uses for S-Ergo Series*

- Ergonomic, lightweight, compact
- Pressure relief, stabilization, and anti-slipage
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider arm pad, providing comfort/support

## S-Ergo 115

### *Standard Lightweight Ergonomic Manual Wheelchair*

- Medicare code K0005
- High strength, weighs only 25 lbs. (w/ footrests) & only 19 lbs in full travel form!
- Anti-Bacterial Upholstery
- Folding backrest
- Fixed armrest, Swing away footrest
- “Tube-in Center” foot-plate, assures better side leg support
- 24” Rear spoke, polyurethane, high profile, flat free wheels
- 16x17 or 18x17 Seat width
- Frame Color: Pearl Silver & Rose Red
- Upholstery: Silver/Black
- Weight Capacity of 250 lbs

items		F/R												Max. User Weight
Width	wheels													

S-Ergo 115	~	16"x17"	7"/24"	16"	17"	23"	12"	8"	39"	20"	16"	36"	25 lbs.	250 lbs.
S-Ergo 115	~	18"x17"	7"/24"	18"	17"	25"	12"	8"	39"	20"	16"	36"	25 lbs.	250 lbs.

*\*Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery*  
*\*\*All models come equipped with 7"x1" front polyurethane wheels*



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue. City of Industry. CA 91748-1418

Tel.626-581-2235 Fax.626-581-2335