



The ultimate in
safe, secure bathing

Bathmaster
sonaris 2

Specifications

Seat Size (including side flaps)	20" x 18" x 2.6"
Seat Size (excluding side flaps)	20" x 15"
Seat Depth	20"
Maximum Seat Height	18"
Minimal Seat Height	2.6"
Baseplate Footprint	21" x 9"
Backrest Size	24" x 14"
Maximum User Weight	375 lbs.
Maximum Recline Angle	40°

0816-21-796—Sonaris2 Bathlift with Seat and Back Covers

Accessories

Swivel Seat Supports the turning movement of getting the legs into the bathtub.

Swivel Transfer Seat Serves to turn and slide the body to the center of the lift.

To order: www.pattersonmedical.com
1-800-323-5547

Great
performance
at a great price



Bathmaster
deltis

Specifications

Seat Size (including side flaps)	19.75" x 28"
Seat Size (excluding side flaps)	19.75" x 14.75"
Seat Depth	20"
Maximum Seat Height	18"
Minimal Seat Height	2.6"
Baseplate Footprint	21" x 9"
Backrest Size	26.75" x 13.75"
Maximum User Weight	309 lbs.
Maximum Recline Angle	40°

0816-64-590—Deltis Bathlift with Seat and Back Covers

Accessories

Swivel Seat Supports the turning movement of getting the legs into the bathtub.

Swivel Transfer Seat Serves to turn and slide the body to the center of the lift.

To order: www.pattersonmedical.com
1-800-323-5547

Bathmaster
sonaris 2



Bathmaster
deltis



Just the lift you need for
easy, independent bathing

Exclusively from Patterson Medical


PATTERSON
MEDICAL

081688167



Bathmaster
sonaris2™

What makes Sonaris2 unique?

Combines Quick Charge convenience, the latest technology, and a higher weight capacity for the ultimate in safe, secure bathing

Super Fast Quick Charging— Advanced battery fully charges in only one hour! 5 minutes is all it takes to charge device for one complete lift (down and up).

Nanotechnology Antimicrobial Protection—Not just a coating, silver ion antimicrobial protection is infused into easy-clean, plastic for long-lasting performance without the spread of bacteria.

Easy Assembly—Seat and back click together or snap apart.

Intelligent Design—Hand controller is triple-sealed for safe submersion and will not lower the unit if battery power isn't sufficient to raise the user out of the tub.

Weight Capacity—375-lbs.

Both Sonaris2 and Deltis feature these impressive details



Bathmaster
deltis™

What makes Deltis unique?

Great performance at a great price and one of the lightest reclining bathlifts on the market today

Super Fast Quick Charging—Advanced battery fully charges in only one hour! 5 minutes is all it takes to charge device for one complete lift (down and up).

Lightweight Design—Bathlift with handset and battery weigh only 24 lbs.

Simple Assembly —Separates easily into two parts for easy installation, removal and transport.

Warning Light— Lets you know when charging is required.

Weight Capacity—309-lbs.



Built-in Safety Controllers float and are fully submersible. Suction cup feet hold bathlifts securely in place.

Compact & Convenient All-plastic design is corrosion-free, easy to clean, and lightweight for effortless transport and storage. Bathlifts fit in virtually any tub and large side flaps make for easy transfers.



Independent bathing has never been easier or more comfortable!

