

**VivaLIFT!**<sup>®</sup>  
POWER RECLINERS  
By Pride<sup>®</sup> Mobility Products Corp.

THE MOST  
COMFORTABLE  
RECLINER  
IN THE WORLD!

## THE 4-MOTOR RECLINER

WITH  
COUNTLESS  
POSITION  
COMBINATIONS



Shown in Stonewash  
Mossy Fabric

# SIERRA

We're always on the go, moving from one activity to the next. Relish those few moments of comfort and closeness when time slows down. The Sierra Collection by **VivaLift!**<sup>®</sup> delivers lasting comfort. Infinite, lay-flat positioning maximizes your comfort while connecting with your loved ones. Independent functions on the Sierra recliner allow you to put your feet up with the footrest extension, while keeping the power headrest in the precise position to support your neck.



**Pride Mobility**<sup>®</sup> (US) 800-800-8586 | (CAN) 888-570-1113  
pridemobility.com | prideliftchairs.com |    

# SIERRA

COLLECTION	SIERRA v.2
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift
WEIGHT CAPACITY	375 lbs.
SEAT-TO-FLOOR	19.5"
SEAT DEPTH	20.5"
SEAT WIDTH	21.5"
TOP OF BACK TO SEAT	25"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	44.5"
DISTANCE FROM WALL	16"
OVERALL WIDTH	35"
CHAIR WEIGHT	144 lbs.
WARRANTY	7-Yr. Pro-rated

## FUNCTIONS AND FEATURES:

-  POWER HEADREST
-  POWER LUMBAR
-  USB REMOTE
-  LITHIUM BATTERY BACKUP
-  FOOTREST EXTENSION
-  POWER BACKREST

## POSITION GUIDE:

-  RELAXED
-  RECLINED
-  UPWARD LIFT
-  INFINITE LAY FLAT

## STANDARD FABRIC: STONEWASH



MOSSY

AMBER

DUNE

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.



**Pride Mobility**® (US) 800-800-8586 | (CAN) 888-570-1113

pridemobility.com | prideliftchairs.com |    