

www.ameriglide.com  
800-987-6308

**HORIZON PLUS**

straight & simple



**AmeriGlide**  
Accessibility  
Solutions

# HORIZON PLUS

straight & simple

The AmeriGlide Horizon straight stairlift is designed for ease-of-use and fast, hassle-free installation. A Horizon stairlift will make difficult staircases accessible again and grant you the freedom of your house once more.

The compact design and floor-mounted straight rail ensure the Horizon takes up as little space as possible and does not impact upon your home's decor.



Easy to use joystick



Secure key operation



Digital diagnostic display



Optional powered hinge



Multiple Seat Colors

- Fixed seat and seat pad
- Manual swivel
- Maximum user weight - 308 lbs.
- Anodised aluminium single rail
- Rail angles 25° to 62°
- Digital diagnostic display
- Infra-red remote controls
- Retractable lap seat belt
- Outdoor model available
- **Optional** Powered hinge rail



Standard Seat



Ergo Seat



## WARRANTY

2 Years: Component parts including the circuit boards

5 Years: Drive train (motor and gear box), gear rack and pinion gear