

STAIRCASE MEASURING GUIDE

for Straight Stairlifts

In order to ensure that your Straight Stairlift installation goes smoothly, it is imperative that proper and accurate measurements are provided. If you have any questions, or would like some assistance with measuring your staircase, please call **Customer Service**.

1. Overall Length

ft. / in.

Stand at the top of your stairs, extend your tape measure until it is laying on all the stair treads and the tip of the tape measure is touching the landing at the bottom of the steps.

2. Top Stair to Bottom Stair

ft. / in.

Stand at the top of your stairs, extend your tape measure until it is laying on all the stair treads and the tip of the tape measure is touching the landing at the bottom of the steps.

3. Staircase Width

ft. / in.

Measure the width of your staircase and enter the measurement.

4. Bottom Clearance

ft. / in.

If there is a wall or a door at the bottom of the steps, measure from the bottom step to the wall.

5. Rise of One Step

ft. / in.

Measure from the top of a step to the very next step down.

6. Tread Depth

ft. / in.

Measure from the front of one step to the front of the kickplate.

7. Diagonal

ft. / in.

Measure from the front of one step to the front of the next step down.

8. Top Clearance

ft. / in.

Measure from the top step to the door, wall, or other obstruction at the top of the stairs.

