




### Standard Fabrics:




Calypso



Copper



Hazelnut



Shiraz



Sterling

This fabric stands out for its comfort, smooth texture, and rich color options. It features our exclusive SofTouch™ with Stain Defense coating.

\* An upcharge will apply for all crossover fabrics.

## EZ Sleeper PR735-MLA



*Comfort Zone*™  
by GOLDEN

The PR735 EZ Sleeper is just for you! This luxury power lift recliner combines a contemporary design with our patented MaxiComfort® positioning for unmatched style and comfort. Its ultramodern, clean lines pair perfectly with a contemporary décor, especially in smaller spaces. Enjoy the long-lasting comfort of this rounded, track arm style chair with a tailored seam back. Your choice of luxurious fabrics and MaxiComfort® positioning will make this chair the favorite seat in the house.

### Standard Features:

- Patented MaxiComfort® Positioning
- Programmable AutoDrive™ Hand Control
- Extra Pocket



# MAXICOMFORT® Series PR735

Model Number	PR735-MLA
Size	Medium/Large
Weight Capacity	375 LB
Positioning	MaxiComfort
Comfort Zone(s)	2
Twilight Technology	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Bolt
Head Rest Cover	No
Arm Covers	No
Overall Width	34.5"
Overall Height	45.5"
Width Between Arms	21.5"
Floor to Top of Seat	20"
Seat Depth	22.5"
Seat to Top of Back	30"
Seat to Top of Arm	7.5"
Distance Required from Wall Reclined	20"
Extension Reclined from Back of Seat to End of Foot Rest	47"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Hand Control Routing	Inside Arm
Extra Pocket Option	Standard
Heat & Massage Option	Yes
Heat Wave Technology Option	No
Articulating Headrest Option	No
Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes
Straight Lift Option	Yes
USB Charging Port	No
Battery Backup	No
Place of Origin	USA



Shown in Shiitake\*

SmartTek Diagnostic System STANDARD on all Golden models.




## Comfort Zone™ by GOLDEN

**Comfort Zones 1 & 2** work together to create individualized wellness positioning and comfort. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive hand control.



**WE PATENTED COMFORT™**



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374  
Fax: 800-628-5165  
www.goldentech.com



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with  
UL 73 and  
CSA C22.2 No. 68



www.facebook.com/goldentechnologies

SHT-PR735-071222