

As Unique as You and Your Home

One of a Kind

The Rave Curved HD Stair Lift is one of a kind. AmeriGlide designs are crafted to fit your needs, staircase, and decor.

The Rave Curved HD Stair Lift is built to last. Every staircase that is fitted with a Rave Curved HD Stair Lift is measured with a calibrated camera kit for digital accuracy. Constructed with highly durable steel, custom precision-bent rails, and finished by experienced craftsmen, this true curved rail is highly engineered to give you an extraordinarily smooth ride.

Your unique path to independence will allow you to enjoy time with your family and your home on every level giving you independence, freedom and the quality of life you deserve.







Tighter Curves

The Rave Curved HD Stair Lift's narrow profile mounts closer to the wall, hugging the stairway with tighter bends and flowing radius curves. This means more room for you, and a comfortable smooth ride.

The Rave Curved HD Stair Lift is Customized for You



20 **Feet Per Minute**

350 **Pound Capacity**

16" **Folded Width**

164' **Max Travel Length**

JANUARY 2021

Rave Curved HD Stair Lift Product Specifications

Weight Capacity Drive System	350 lb (158.8 kg) Steel Rack and Pinion	Warranty	Limited Lifetime Warranty, 2-Year on Parts, 1-Year on Battery
Lift Speed	20 fpm	Options	Emergency Safety, Wide Seat 21" (53.3 cm), Key Switch, Additional Call/Send, Additional Charging Station, Custom Colors
Seat	19" W x 15" D x 20.5" H		
Seat	(48.3 x 38.1 x 52.1 cm)		Charging Station, Custom Colors
Folded Width	16" (40.6 cm)		
Control Circuit	24V DC		AmeriGlide
Standard Color	Beige		AIICIUIUCE ACCESSIBILITY SOLUTIONS
Max Rail Length	164' (50 m)		IANUARY 2021