

WITH P

12.5" (Folded from wall)

I love my Crest
HD Stairlift for
the freedom &
peace-of-mind
it gives me!
?

CREST HT

LIMITED LIFETIME WARRANTY















## Features:

- Folds to within 12.5" of the wall
- Attaches to steps, not the wall
- Narrow, low profile track (2.5" tall x 6.5" wide)
- Footrest obstruction sensor with audio beep notification
- Battery power allows up to 20 trips during power outages
- Full length continuous charge system with smart charge technology
- Rocker switch conveniently located on either armrest
- Constant pressure controls
- Wireless remote controls for upper & lower landings
- Chassis mounted installation load switch
- Quickly installs on left or right side of stairway
- Rack and pinion drive with concealed gear rack
- High seat back
- Low-to-the-floor footrest
- 3" to 6" horizontal extension at top landing
- 25 to 45 degree angle applications
- ETL listed ASME A17.5 / B44 and A18.1
- 350 lbs capacity
- Overspeed governor
- Proudly made in the USA!

85 degree swivel seat at top landing

## **Technical Information**







Dimensions subject to change without notice. Staying Home Corporation • PO Box 37, Harrisonville, MO 64701 info@stayinghome.com • 877-378-4275



Best Warranty In The Industry!

1 Year	warranty on BATTERIES
2 Year	warranty on PARTS
5 Year	warranty on DRIVE TRAIN
	& CIRCUIT BOARD